



 **ECHOMEDIC**

**EB ONE V1**

### **Draagbaar echosysteem Rund- Paard**

Gebruiksvriendelijk !

Mode : B, B/B, B/M, M,A

Transducers : endorectaal linéair , endorectaal semi-convex ,  
microconvex, convex , linéair

Metingen : distance, area ,OB, backfat , lean rate  
Obstetrische meetsoftware rund, paard, etc

Gewicht 1,1kg

Afmetingen : 158X158X65 ( mm, LxWxH)

+ 5 presets

+10 verschillende rooster types

**Autonomie :5 uur**

**5,8GHz RFID** connection 'wireless external monitor'

**Wifi : 802.11g/20Mhz/2.4G**

APP support Android, Windows

Vector : PAD/Phone/PC

**Waterdicht !**

**1 jaar garantie met ECHOMEDIC 7/7 service**

**Vervangtoestel binnen de 48h**

**Eigen technische dienst**

**ECHOMEDIC**

**info@echomedic.be T + 32(0)468/115967**