

**mindray**

## **Z60 Vet**

Diagnostic Ultrasound System

# Wise Choice for Professionals



# Z60 Vet

The Z60 Vet is a powerful and versatile Color Doppler system which provides you the best solution for veterinary ultrasound imaging with excellent performance. Due to the comprehensive configurations, dedicated veterinary workflow and integrated design, you can be assured of an outstanding ultrasound experience in both large and companion animals. Packed with professional innovative technology inside this smart portable system makes Z60 Vet a perfect veterinary ultrasound system with excellent image quality and various advanced features.

## Comprehensive Veterinary Applications

- **Veterinary scanning modes:** B, Color, PW, CW, Power, M, Free Xros M
- **Veterinary transducer variety:** Micro-Convex, Linear, Convex, Phased array, Endocavity
- **Veterinary applications:** reproductive organs, abdomen, cardiac, tendinous, small parts and more
- **Veterinary exam packages:** Measurement, Comment, Report and Body mark designed for canines, felines, equines, bovines, ovines and more



## Powerful Cardivascular Solutions

- **PW/CW Doppler:** Auto Doppler trace and calculation
- **Color M:** includes Color Flow M mode and Color Tissue M mode
- **Auto IMT:** Automatic measurement for veterinary Intima-Media Thickness
- **Tissue Doppler Imaging:** Full functional (TVI, TEI, TVM, TVD)
- **ECG module** for veterinary cardiology analysis




## Dedicated Veterinary Workflow


- **iTouch™:** Quick image optimization in one touch
- **iZoom™:** Instant full screen image expansion with one button
- **User-defined keys:** programmable personalized operations tailored for veterinarian

## Ergonomic Mobile Design

- 60 degree tilting angle adjustable monitor
- Integrated design with internal AC power adapter
- 1.5h uninterrupted scanning with rechargeable battery
- 15 inch High Definition monitor with full screen design
- Water-proof protection cover designed for veterinary keyboard
- Backlit control panel designed for veterinary diagnosis in clinics
- Can be packed in a convenient hand-carried bag for easy transportation



Canine Kidney



Canine Intestine

[www.mindray.com](http://www.mindray.com)

P/N:ENG-Z60 Vet-210285x2P-20190909

©2019 Shenzhen Mindray Bio-Medical Electronics Co.,Ltd. All rights reserved.

**mindray**  
healthcare within reach